

# Mobile Phase Preparation, Storage, & Handling

---






## About Sputnik V

- Developer of Russia's Sputnik V Covid-19 vaccine claims there were no cases of blood clots encountered during trials & use of vaccine
  - ...pointed to a recent [study](#), published in The New England Journal of Medicine, that linked blood clots in those inoculated to insufficient purification of vaccines.
  - *Unlike other vaccines,* the Russian [dose] goes through a four-stage purification process, **consisting of two stages of chromatography** and two stages of tangential flow filtration, which assures its “*quality and safety*,” the statement read.

# Chromatography in Pharma



## GAS CHROMATOGRAPHY

Any time a lab technician must separate volatile materials or identify raw materials in a mixture, many pharmaceutical professionals reach for gas chromatography (GC). Since GC machines are capable of analyzing extremely small and light compounds, GC is used in post-production.

## HIGH-PERFORMANCE LIQUID CHROMATOGRAPHY

Efficient and accurate results are what make HPLC the primary chromatography method in the pharmaceutical industry. HPLC uses a liquid as the mobile phase to ensure the fastest process. Furthermore, HPLC machines have a column as the stationary phase, which separates all the compounds accordingly.



## DISCOVERING NEW MEDICINE

These two chromatography methods have allowed professionals to identify a molecule with remedial properties. After finding a particular molecule—or molecules—lab technicians can begin research by developing new formulas. From there, there will be further research and analysis to “fine-tune” the medicine before launching it to the market.




## ANALYSIS AND IMPROVING AN EXISTING MEDICINE

Pharmaceutical lab technicians are always analyzing existing medicines. By doing so, they can safely alter the formulations to improve results and potentially reduce or eliminate side effects.


## SOME PHARMACEUTICAL COMPANIES WILL COMBINE METHODS

Since each chromatography method has its own strengths and weaknesses, using an additional method can reduce any inconsistencies in findings. In other words, adding an additional method of chromatography into the research process will certainly help formulate a stronger hypothesis. Furthermore, additional analytical instruments provide an additional layer of data to work with.



GENTECHSCIENTIFIC.COM


# Impurity Analysis of Ziprasidone



B.


	Name	Component Type	Retention Time (min)	Area (μV*sec)	Impurity Response	Above Reporting Threshold	Above Identification Threshold	Above Qualification Threshold	ICH Threshold	Pass?
1	Compound A	Specified Impurity	0.724	439	0.027	✓	✓	✓	Above Qualification Threshold	Fail
2	Compound B	Specified Impurity	1.824	230	0.017	✗	✗	✗	Below Reporting Threshold	Pass
3	Ziprasidone	Main Component	1.911	1145524						
4	Compound D	Low Level Impurity	2.978							
5	Compound C	Low Level Impurity	3.119							
6	Sum Specified Impurities	Specified Impurity		669	0.044					Fail

# High Performance Liquid Chromatography




# Creating Separations

- Typically want to have analyte retained somewhat by column
  - Choose a SP with a similar polarity
  - Choose a MP with a dissimilar polarity



# Normal Phase Chromatography

- SP = polar (more than MP)
  - Ex. Silica, cyano, diol amino
- MP = non-polar (less polar than SP)
  - Ex. Hexane
- Polar compounds retain most strongly
- Non-polar elute more easily
- Called “normal” bc it was the first type of partition separation developed



# Reverse Phase Chromatography


- Most common form of liquid chromatography

- HPLC typically refers to RP




- SP = non-polar (more than MP)

- Ex. C18, C8, C4, phenyl



- MP = polar (more than SP)

- Ex. Water, MeOH, ACN




- Non-polar compounds retain most strongly
- Polar elute more easily
- Called “reverse” bc it is the opposite of normal phase





# Demonstration 25-1 Normal- and Reversed-Phase Chromatography (3 of 3)

- The normal-phase separation uses a  $1 \times 6.5$ -cm silica thin-layer chromatography plate and is developed in a  $1.5 \times 10$ -cm test tube using  $\leq 1$  mL of 70:28:2 (vol:vol:vol) heptane:acetone:anhydrous ethanol.
- The reversed-phase separation uses a  $1 \times 6.5$ -cm  $C_{18}$  thin-layer chromatography plate and is developed in a  $1.5 \times 10$ -cm test tube using  $\leq 1$  mL of 15:35:50 (vol:vol:vol) heptane:acetonitrile:95% ethanol.
- With a pencil, gently draw a line  $\sim 0.7$  cm from the end of the plate.
- Use a capillary tube to apply a series of spots along the horizontal line.
- Add eluent to the test tube and place the plate in the test tube (make sure the eluent does not rise above the horizontal line).
- Remove the plate when the solvent front is 0.5–1.0 cm from the top ( $< 10$  min).



(a) Normal-phase (b) Reversed-phase  
Courtesy of L. Kvittingen, Norwegian University of Science and Technology. From B. J. Sjørnes, L. Kvittingen, and R. Schmid, "Normal and Reversed-Phase Thin Layer Chromatography of Green Leaf Extracts," *J. Chem. Ed.* **2015**, 92, 193. Reprinted with permission © 2015, American Chemical Society.

# Simplified RP Retention Mechanism

- Retention based on hydrophobic interactions
- The high polarity of the water “repels” the non-polar analytes out of the MP and onto the SP
  - The concentration of the analyte on the SP increases





Figure 1: Schematic representation of Reversed Phase HPLC



# Determination of Elution Order

- Hydrophilic/polar compounds will elute first




# Elution Order

---

- The structure of a sample gives clues as to elution order
- When comparing SIMILAR compounds
  - Less water soluble = more retention
  - More carbon atoms = more retention
  - Branched chain compounds = less retention
  - Unsaturation = less retention
  - Polar and charged species = less retention
- These are GENERALIZATIONS
  - SP and MP dependent

# General Trends in Elution






Note this is a homologous series!




# Practice

- Indicate the elution order by RP (1=first)
- Butyl, ethyl, methyl, propyl paraben



# High Performance Liquid Chromatography



# Solvent Strength


---

- When we make the mobile phase more similar to the stationary phase we increase its strength
  - NP – make it more polar
  - RP – make it more non-polar
- The stronger the solvent, the faster compounds elute
  - Regardless of analyte properly




# Mobile Phases in RP

- Water is least like the SP
  - Typically termed the “weak solvent”
- MP that are more like the SP are “stronger”
  - Elution faster




# Strength of the Mobile Phase

- We can change the “strength” of the MP in two ways for RP separations
  1. Increasing the % organic in our mixture
    - 30:70 MeOH:water vs 70:30 MeOH:Water
  2. Changing to a less polar solvent
    - Water to MeOH or ACN or THF




# RP Isocratic – Increasing % Organic



Not all of the peaks are resolved

Peak identification

1. Pindolol
2. Diisopyridamide
3. Propranolol
4. Dipyridamole
5. Diltiazem




- Good resolution
- Fast analysis
- No time wasted




Increasing % organic reduces retention time

- Organic is the strong solvent







# Change in Solvent Strength



- B = organic, RP separation
- Increasing organic reduces retention time – stronger solvent

# Change in Solvent Strength

B=organic, RP separation



# Impact of Increased % Organic

- Why reduced retention time?
  - More organic makes the MP more like the SP
  - Less drive of the analyte onto the SP




Figure 1: Schematic representation of Reversed Phase HPLC



# Impact of Increased % Organic (RP)

---

- **ALL** compounds elute earlier
- Increase organic by 10% can decrease  $k$  by a factor of 2-3
- Typically a change of 5 to 15 % in retention time per 1 % change in organic composition
- This is a problem for reproducibility






# HPLC Simulator

---


- [https://www.multidlc.org/hplcsim/4\\_2\\_0/](https://www.multidlc.org/hplcsim/4_2_0/)
  - Look at the compounds benzylformamide, benzylalcohol, phenol, acetanilide, and phenylpropanol
  - Use water and ACN as the mobile phase with an isocratic method.
  - Change from 20%B to 22% B?
    - What happens?
  - Take a few minutes and play with the mobile phase conditions, the gradient options and the column properties - observe what happens to the chromatogram as you change these parameters



# Impact of More Non-Polar MP



- Identical conditions, just ACN in lieu of MeOH
- ACN is less polar than MeOH
- ACN = stronger solvent
  - Moves things off the column faster



# Impact of Stronger Solvent

- Why reduced retention time?
  - More non-polar organic makes the MP more like the SP
  - Less drive of the analyte onto the SP




Figure 1: Schematic representation of Reversed Phase HPLC




# Organic Modifier

---


- Methanol

- Protic solvent – hydrogen bond donor
- Weaker elution solvent
- Higher viscosity than acetonitrile




- Acetonitrile


- Aprotic solvent – hydrogen bond acceptor
- Strong elution solvent
- Low viscosity




# Changing Organic Modifier




- Changing solvent or “organic modifier” can also change selectivity
- Must always confirm the elution order if make a change



# RP Isocratic – ACN to MeOH



- MeOH is more polar than ACN → it is a weaker solvent
- Weaker = less able to compete with the stationary phase
  - Greater retention by reverse phase



- May also be some selectivity differences




# Impact of Organic

---

- Can use this knowledge to
  - Develop or modify methods
    - Elute faster or differently
  - Troubleshoot unexpected results

# Troubleshooting

- You are working in A240 running the parabens experiment, your method calls for 70:30 MeOH:water and you expect this elution order and time. Instead you find
  - The elution time for all peaks is later than expected, what should you check?



# Troubleshooting


---

- You are running a published method using 50:50 water:MeOH
- The method worked yesterday (!) but now you find that most of your peaks are eluting earlier AND two have swapped places in the elution order
- What should you check/try?




# Protonation State


- Predict the protonation states of the following at pH = 4, 7, and 10




Phenol  
pKa = 10




Aniline  
pKa = 10




Triethylamine  
pKa = 10.75




amphetamine  
pKa = 9.9




Alanine  
pKa-COOH = 2.34  
pKa-NH3 = 9.69



ASA  
pKa = 3.49




Lysine  
pKa COOH = 2.16  
pKa-NH3 = 9.06  
pKa - side chain NH3 = 10.54



Butanoic  
acid  
pKa = 4.82



# Impact of Mobile Phase pH



# Impact of pH on Retention





- Compounds:

- p-Toluamide – neutral
- Lidocaine – base, pKa 7.8
- Ibuprofen – acid, pKa 5.2




# Mobile Phase: Buffers

- When would we use a buffered mobile phase?
  - When pH is important
  - What types of molecules would be involved – Ionizable or un-ionizable?
    - Ionizable





# Effect of pH on Retention

## Reversed Phase Retention Map



# Effect of pH on Retention



Silica-based materials dissolve at high pH



# Consistency and pH

---

- Need to control pH to ensure consistent results
  - Especially if near the pKa of the analyte
- Ideally want to be 2 pH units from the pKa

# Are Buffers Always Used?



---

- No
- In some cases Trifluoroacetic acid (TFA) or formic acid (FA) are added in small amounts (0.1% v/v) to the mobile phase
  - Pull the pH down, assumed that exact pH not a factor



# Reverse Phase and pH

- Some molecules can have multiple charged states



- Reverse phase methods use a non-polar surface
- What is will retain better by RP?
  - Charged or uncharged?
- How can we keep it in one state or another?
  - Through controlling the pH of the mobile phase




# Troubleshooting

---

- Your solution to be analyzed contains
  - An amine ( $pK_a = 9$ ) and other un-ionizable species
- Your conditions
  - 50:50 MeOH:pH 11 buffer (aq)
- You run your analysis and have excellent separation
- The next day you make up new mobile phases and find that the amine has reduced retention but all other peaks are unchanged.
  - What might be the problem?

# High Performance Liquid Chromatography



# USP: Liquid Chromatography

---

From General Chapter <621>

The term "liquid chromatography" (LC), as used in the compendia, is synonymous with high-pressure liquid chromatography and high-performance liquid chromatography.

LC is a separation technique based on a solid stationary phase and a liquid mobile phase.

# USP: Mobile Phase

---

From General Chapter <621>

The mobile phase is a solvent or a mixture of solvents, as defined in the individual monograph.

# USP – Mobile Phase Examples

---

Mobile phase can be simple:

## *Acetaminophen Tablets:*

Assay is HPLC based, Mobile phase is

- Solution A: 1% v/v glacial acetic acid in water
- Solution B: Methanol

How the mobile phase is prepared and stored impacts the results

Mobile phase can be more complicated:

## *Atorvastatin*

Assay is HPLC based, Mobile phase is

- Buffer: 3.9 g/L of ammonium acetate in water. Adjust with glacial acetic acid to a pH of  $5.0 \pm 0.1$ .
- Solution A: Acetonitrile, stabilizer-free tetrahydrofuran, and *Buffer* (21:12:67)
- Solution B: Acetonitrile, stabilizer-free tetrahydrofuran, and *Buffer* (61:12:27)

# The Mobile Phase

---

- Issues with mobile phase cause up to 70% of all LC problems
- What might go wrong with mobile phase?
  - Composition
  - Particulates
  - pH
  - Contaminates
  - Bacterial growth
- We must take great care in preparing, and documenting the preparation of, our mobile phase



# Glassware

---

## Choose your glassware wisely

- Only make up with clean glassware that has not been washed with detergent
  - Detergent can stick to glass and then leech off into the mobile phase, contaminating it
  - RINSE glassware with HPLC grade solvents and water only



<https://mygoldfishisalive.com/why-do-my-goldfish-keep-dying>






# Impact of Detergent

---

- The detergent can then elute off the column during your gradient run causing an unknown or “ghost” peak



- What might be the implications if a ghost peak is seen in during a QC run for batch release testing?
  - Investigations, release delays, \$\$\$\$\$\$\$\$



# Solvent and Water Quality

---

- Choose the highest grade of solvent possible
  - HPLC Grade for HPLC work
  - Optima or similar grade for UPLC and/or MS work
- Choose the highest grade of water possible – Industry standard is ultrapure
  - Undergone successive steps of filtration and deionization & characterized as having a resistivity of  $18.2 \text{ M}\Omega\cdot\text{cm}$  at  $25^\circ\text{C}$
  - Commonly called “Milli-Q” water
    - Term a trademark of Millipore
    - Similar to Kleenex for tissue



# Filtration


---

- Everything that will enter a chromatographic system must be filtered
  - By you or by the manufacturer
  - This includes solvents, mobile phases and standards and samples in solution
- Filtration is happening continuously as the system is being used



# In-Line Filtration


- Performed on solvents and samples that have entered the system, prior to hitting the column
- Filters are placed
  - In the solvent, at the end of the solvent lines
  - In the pump, before the check valves
  - Between the pump and the injector
  - The injector and the column



# External Filtration

---


- On samples prior to injection
- Pore size of membrane filters
  - 0.45  $\mu\text{m}$  pore size for 3+  $\mu\text{m}$  particle size columns
  - Particle sizes <3  $\mu\text{m}$  use 0.22  $\mu\text{m}$  pore size



# Impact of Filtration


## Results – Filtration impact on sub-2 micron column A by Latex Bead 0.3 µm solution

Effects on Filtration on sub-2 micron Column Life



## Results – Filtration impact on sub-2 micron column B by Human Plasma PPT Extract

Effects of filtration on sub-2 column life time



# Filter selection


---

- Filters are available in a variety of physical sizes and membrane types
- Filters must be chosen with their intended use/application in mind and must be used properly



# Membrane Selection

## Agilent Captiva Syringe Filter Selection Guide



- Defined by sample and solvent
- Choosing the wrong filter type can cause issues with
  - Solvent compatibility
  - Sample contamination
  - Sample loss





# External Filtration - Solvents

---

- Based upon application and mobile phase components
- Milli-Q water can be used directly
- HPLC and higher grades of solvents can typically be used from the bottle



# External Filtration - Buffers

---

- Any time a solid is added to the solvent it must be filtered before being put on the system
- Remember: Buffers are used to control the pH of a solution
- All additives should be of the highest grade possible
- Must filter to remove particulates



# Buffers - Salts

---

- Different buffers will control pH changes at different pH ranges
  - Related to their pKa values
- In older methods Na & K salts were very common as were phosphate buffers
  - BUT these are not amenable to MS
- Modern methods tend to move away from non-volatile buffers
  - Ammonium with acetate, formate and bicarbonate are now more common



# Volatile Buffers

---

- Non-volatile buffers interfere with the performance of mass spectrometers
- LC/MS is becoming more and more prominent
  - Environmental – emerging contaminants
  - Biotech/biopharm – important for characterization
  - Pharma – impurity identification
  - Food – residue analysis
  - Useful for compounds without a chromophore
- Often labs without a MS develop amenable methods just in case...



# Drawbacks of Filtration

- For **trace analyses** involving sensitive MS detectors filtration is often avoided
  - Can contribute contamination greater than analyte concentrations
- Often solid additives are avoided for MS work
  - Use high quality liquid additives

Diquat/Paraquat Recovery Data (n=7) for Three Types of Water Spiked at 40 ng/L

Sample Type	Recovery (%RSD)	Recovery (%RSD)
	Diquat	Paraquat
Groundwater	82 (8)	89 (16)
Tapwater	84 (5)	88 (8)
River water	83 (4)	89 (19)




Figure 4. Typical LC/MS(MS) Calibration Curves.



# Mobile Phase Preparation

---

- If we are making a mixed mobile phase (1:1 ACN: 10 mM sodium acetate, pH 4), when do we measure pH?
  - Before we add organic or after?
  - Remember: The presence of organic solvent in a solution affects the accuracy of a pH meter
  - The addition of organic solvents also changes the pH of a solution if not buffered
- Depending on the analytes, minor changes in pH can impact results
- Most modern methods mix solvent ON the system
  - Using line A and B



# Kahoot

---

- Mobile Phase part 1

# Mobile Phase - Degassing

---

- Solvents equilibrate with gases in the atmosphere and the solubility of gases in a pure solvent is higher than that of a mixture
  - When solvents are mixed together, the solubility decreases, the solutions are over saturated and gas bubbles form




<https://www.forbes.com/sites/kalevleetaru/2017/12/18/why-was-2017-the-year-of-the-filter-bubble/?sh=319c4ebd746b>





# Mobile Phase - Degassing

- When do solvents mix in HPLC?



# Mobile Phase - Degassing


---

- Solvents equilibrate with gases in the atmosphere and the solubility of gases in a pure solvent is higher than that of a mixture
  - When solvents are mixed together, the solubility decreases, the solutions are over saturated and gas bubbles form
  - Gas bubbles can occur on mixing or anywhere within the system where a nucleation site exists



<https://www.forbes.com/sites/kalevleetaru/2017/12/18/why-was-2017-the-year-of-the-filter-bubble/?sh=319c4ebd746b>

## Problems with Formation of Bubbles in Flow Lines



## Detrimental Effects of Dissolved Air (Oxygen) on Detection



# Mobile Phase – Degassing Benefits

---

## 1. Baseline stability

- Detector is free of bubbles – no sudden changes in signal & pump is consistent in solvent delivery

## 2. Injection volume reproducibility

- No air in lines to affect the introduction of the band of sample in to the MP stream


## 3. Retention time reproducibility

- Pump operation, solvent composition and sample introduction are consistent run to run



# Mobile Phase – Degassing Methods

- Modern LC systems have built in degassers which are effective for most applications



<http://www.chromatographyonline.com/mobile-phase-degassing-what-why-and-how>



# Mobile Phase – Mixed Solution

---

- Some methods will require solvent mixtures as the mobile phases
  - Particularly older methods
- Examples
  - A = 95% ACN, 5 % water; B = 1:1 ACN:water
  - A = water, B = 90% MeOH, 10% THF
- How the solution is made is critical for success




# Mobile Phase – Mixed Solvent Prep

---

- How would you prepare 1L of a 60:40 MeOH:Water solution?
  1. Measure out 600 mL of MeOH in a graduated cylinder and fill to 1L with water
  2. Measure out 400 mL of water in a graduated cylinder and fill to 1L with MeOH
  3. Measure each separately and combine
  4. Weigh each out separately and combine



# Mobile Phase – 60:40 Mix Prep




Options 1 and 2: Adding organic/aqueous directly to measuring vessel  
Need to specify which solvent is added first: e.g.  $\text{H}_2\text{O}$  to MeOH or MeOH to  $\text{H}_2\text{O}$





# Mobile Phase – 60:40 Mix Prep



**Options 1 and 2: Adding organic/aqueous directly to measuring vessel**

**Need to specify which solvent is added first: e.g. MeOH to H<sub>2</sub>O or H<sub>2</sub>O to MeOH**

**Option 3: Separate volumetric measurement \*\*preferred method\*\***



# Mobile Phase – 60:40 Mix Prep


---

- Measuring out the solvents separately THEN mixing is the preferred method
  - Volumetrically is most common
  - Weighing is most accurate
    - Eliminates any volume differences due to temperature
    - Rarely done
- You must clearly state how your mobile phase is prepared





# Impact of Improper Mixture

- In reverse phase, what is the impact of too much organic? (ACN or MeOH)
  - Things elute too fast



# Remember Strong vs Weak Solvents

- Strong solvents are more like the stationary phase
  - Push analytes off faster
  - Better able to compete for analytes
- Weak solvents are different from the stationary phase
  - All analytes will move slower



# Mobile Phase -Storage

---

- Mobile phase should be covered and insulated to help ensure that its intended characteristics are maintained
- What might change?
  - Composition
  - Contaminants
- MP's containing buffers and water have a finite life.
  - Should be discarded daily or every other day
  - Why?
    - Grow microbes





# Mobile Phase - Storage

---


- Aqueous solutions should be made every day or every other day
  - Why?
  - Aqueous solutions will grow algae and bacteria which will enter the system and clog lines and columns
    - Over pressure and/or peak shape problems

Scanning Electron  
Microscopy (SEM)  
Images of  
Column Frits




# Mobile Phase – Microbial Growth

- Microbial growth causes
  - Blockages in the filters, frits, columns and pump
  - Ghost peaks
    - Especially during gradient elution
  - Increased background noise
- Mobile phases should always be labeled and dated
  - When made and expiration date




Scanning Electron  
Microscopy (SEM)  
Images of  
Column Frits



# Microbes in Mobile Phase

---

- To prevent microbial growth in an idle HPLC system
  - Flush buffers out with water
  - Store the system in minimum of 10% organic (MeOH or ACN)
- Buffers must be flushed from system after use
  - May precipitate and create significant blockage as well as degradation of the piston seals (seals fall apart)





# Mobile Phase - Storage

---

- Ensure sealed to avoid evaporative loss
  - Affects composition
  - Make up fresh every day or every other day
- Do NOT use Parafilm
  - Dissolves in organic solvents
  - Bleeds into solvent
  - Builds up on the head of the column and elutes as a blob-o-gram
- All mobile phase should be labelled
  - When made, by who, what it is, when it expires




# Kahoot

---




# Practice Problems



(a)

- You are provided the data to the left.
- If you were to develop a method, what pH would you choose to have retention of all three compounds with robustness?
  - Why?
- What would you add to the mobile phase to help with robustness?
- What key things would you need to do to ensure uptime for the instrument when preparing and storing the mobile phase?


# Practice Problems



- Explain the difference between pH 4.4 and pH 3.

# Practice Problems

- You have been given a method to validate which separates 4-methylaniline and codeine
  - Column – Waters XBridge C18, 4.6 x 100 mm, flow rate 1 mL/min
  - MP = 50:50 MeOH:10 mM sodium bicarbonate, pH 8
- Do you have any concerns do you have about the robustness of the method
  - Why?



(b)

# Summary

---

- The preparation of our mobile phase is key for chromatographic success
  - Glassware – detergent?
  - Water and solvent quality
  - Filtration – compatibility and porosity
  - Proper mixing – pH before, measure each component separately
  - Degassing – essential for good performance
- Proper storage is also essential to avoid evaporative loss and bacterial growth
- Know the implications of each of these

