

Troubleshooting & Chromatographic Terms

Troubleshooting

- When you run into a problem, what are good steps to follow to solve the problem?
 - Ex: Trying to win a video game?
 - Ex: Trying to fix a printer issue?
 - Ex: Trying to fix a broken car?

Troubleshooting 101

- In Chromatography problems can arise from many sources
- Effective troubleshooting means
 1. Defining the problem
 2. Isolating the source (finding the cause)
- How should effective troubleshooting be done?
 - Logically, systematically, intelligently, experientially, collectively
 - This is true for all troubleshooting

Troubleshooting 101

1. Don't add to the problem

- Think before you act
- Once the problem is evident, make sure that your actions will not contribute to the problem: don't make things worse

2. Define the Problem

- Identify the symptoms and accurately define the problem
- The real cause may not be the most apparent thing
- Only by identifying the true problem can you find the correct solution
- Symptoms are not the same as problems
- **SIMPLE THINGS FIRST**

Troubleshooting – Define the Problem

- What is the problem/symptom?
 - Runny nose
 - Cough
 - Aching joints
- These are all symptoms, the cause is the flu

Troubleshooting Instrumentation


Need to know how each piece fits together to properly troubleshoot the instrument

The HPLC system consists of

- Solvent
- A solvent delivery system (pump)
- A sample injection valve
- A chromatographic column
- A detector
- A data system


Other important aspects to consider

- Sample and standard preparation



Chromatograms as Tools

- Used to identify and quantify



Chromatograms as Tools

- Can also tell us about
 - The health of the column and the instrument
 - The suitability of the instrument & column for a particular analysis





Image credit: ChromAcademy

Chromatography Terms

- We can describe a separation using a number of parameters including
 - The separating power of a column
 - The retention of a compound on a column
 - The degree of separation between two peaks
- These terms allow us to compare different analyses, columns, instruments, and methods

Retention Time


- Chromatogram is a graph of detector response (y) vs time (x)
- Retention time t_R
 - How long it takes for a compound to go from the injector to the detector



Harris/Lucy, *Quantitative Chemical Analysis*, 10e, © 2020 W. H. Freeman and Company

Mobile phase or an unretained solute (methane) travels through the column in the minimum possible time, t_M .

Void Time/Void Marker




Void time t_0 or t_M

- The time it takes an **unretained species** to reach the detector after being injected
- Void markers are included in a sample or injected separately to determine void time



Retention Time: t_R

- Chromatogram is a graph of detector response (y) vs time (x)
 - Retention time t_R
 - How long it takes for a compound to go from the injector to the detector
 - Elution time
 - Adjusted retention time t'_R
 - The time it takes for a compound to elute, beyond that required for the void marker
- $$t'_R = t_R - t_M$$



Harris/Lucy, *Quantitative Chemical Analysis*, 10e, © 2020 W. H. Freeman and Company


Mobile phase or an unretained solute (methane) travels through the column in the minimum possible time, t_M .

Retention Factor: k

- Retention factor
 - Used to describe the degree of retention for a compound

$$k = \frac{t_R - t_M}{t_M} = \frac{t'_R}{t_M}$$

- As $t_R \uparrow$, $k \uparrow$



Harris/Lucy, *Quantitative Chemical Analysis*, 10e, © 2020 W. H. Freeman and Company


Selectivity or Separation Factor: α

- α is the ratio of the adjusted retention times for two solutes.
- Provides relative selectivity of compounds
 - Often used as a criteria for a method

$$\alpha = \frac{t'_{R2}}{t'_{R1}} = \frac{k_2}{k_1}$$

where component 2 is eluted later than component 1

Figure 23-7



Harris/Lucy, *Quantitative Chemical Analysis*, 10e, © 2020 W. H. Freeman and Company

Practice Problem 2

- Compounds A and B have retention times of 2.142 and 2.193 min, respectively, on a 15-cm column.
- An unretained species passes through the column in 0.864 min.
- The peak widths for A and B are 0.051 and 0.052 min, respectively.
- Calculate:
 - The separation factor/selectivity for the separation

$$\alpha = \frac{t'_{R2}}{t'_{R1}} = \frac{k_2}{k_1}$$

Selectivity or Separation Factor: α


- The term selectivity is also used qualitatively to mean the elution order and relative degree of separation for a set of compounds

Selectivity Changes - Example

Conditions :

Columns: 2.1 x 100 mm


- We would say the selectivity is changes for these different columns
 - Meaning the elution order and/or relative degree of separation



- 1,3,5-TNB
- 1,3-DNB
- NB
- Tetryl
- TNT
- 2-Am-4,6-DNT
- 4-Am-2,6 DNT
- 2,4-DNT
- 2,6-DNT
- 2-NT
- 4-NT
- 3-NT


Efficiency or Plate Count

- N = the efficiency of a column
 - A gauge of how well the column is performing
- Should stay constant
 - If changes perhaps the column is damaged or not connected properly
- $W_{0.5}$ is easiest to measure
 - Esp if baseline separation not achieved





Practice Problem 3

- Compounds A and B have retention times of 2.142 and 2.193 min, respectively, on a 15-cm column.
- An unretained species passes through the column in 0.864 min.
- The peak widths at base for A and B are 0.051 and 0.052 min, respectively.
- Calculate:
 - The efficiency for the last peak



Resolution Between Peaks

- Resolution (R_s) = $\frac{2 * |t_{RB} - t_{RA}|}{W_A + W_B}$
- Separation at **base** of peaks
- Baseline resolution is ideal
 - Where the peaks touch down between each other ($R_s = 1.5$)
- Many methods have a minimum resolution requirements



Practice Problem 4

- Compounds A and B have retention times of 2.142 and 2.193 min, respectively, on a 15-cm column.
- An unretained species passes through the column in 0.864 min.
- The peak widths for A and B are 0.051 and 0.052 min, respectively.
- Calculate:

- $Resolution (Rs) = \frac{2 * |t_{RB} - t_{RA}|}{W_A + W_B}$


Selectivity vs Resolution

- Both take into account the separation of the peaks at their apex
 - Only resolution includes the width of the peaks


Tailing Factor

- In an ideal world peaks are symmetrical
- In HPLC “tailing” can be problematic
 - Decreases resolution
 - USP tailing calculated by software
 - Many methods have a USP Tailing factor pass/fail criteria



USP Tailing Factor



Asymmetry




Examples of Tailing



Applications to Troubleshooting

- We can compare these terms to monitor for changes and evaluate methods for success
 - Changing retention times may indicate problems with the flow rate, column, or mobile phase
 - Changing resolution may indicate problems with tubing and/or columns
 - Changing peak areas may indicate a problem with the injector, samples, or standards
 - Changing efficiency may indicate the column is aging and in need of replacement
 - Unacceptable values may indicate the instrument is not suitable or that the method is not fit for purpose

Troubleshooting - Chromatograms



“My peak areas aren’t the same”


- Define the Problem
 - All of the peaks or just some?
 - How many injections?
 - How significant?
- *All peak areas for second injection were double the first*

Troubleshooting 101

3. **Gather information**

- What is and is not affected?
- Do I know and understand these things?
- What was happening before the problem happened?
 - What events transpired
- Give information that is accurate, detailed and all inclusive
- Don't leave something out because you think it does not matter
- Review critically and carefully any hard data the supports the existence of the problem

Troubleshooting - Chromatograms




- *Peak areas for second injection were double the first*
- Gather information
 - What is the standard? Or sample?
 - Same instrument? Manual or Autosampler?
 - Same day? Right after each other?
 - Retention times unchanged
 - No new peaks

Troubleshooting 101

4. **Narrow the cause down**

- Rule #1 – consider the simple things first
- Eliminate potential causes that would not apply
- To investigate the issue change or consider ONE thing at a time

Troubleshooting – Chromatograms




- Potential causes
 - User error with needle (manual injector)
 - Variability in needle fill
 - Needle not fully seated in injector
 - Injector malfunction (autosampler)
 - Wrong standard injected

Troubleshooting 101

5. Try the options

- What is left to try & does it make sense to try it?
- Keep the problem and what is and is not affected in the front of your mind
- Get help
 - Choose your help wisely!!!

Troubleshooting – Chromatograms




- Information: Autosampler, making multiple injections in series
- Options: Wrong standard injected, autosampler malfunction
- Try: Review injection instructions – did you inject from the same vial twice?
 - Oh! The method says vial 1 then vial 2!

Troubleshooting 101

6. **Evaluate success**

- Did your solution work?
- Go tell your boss what happened and how you took care of it:
Congratulations - you are gaining experience and prestige!

Troubleshooting – Chromatograms



- Updated method to inject twice from the vial 1 before moving to vial 2.
 - Problem solved

Troubleshooting 101 - Summary

1. Don't add to the problem

2. Define the Problem

3. Gather information

4. Narrow it down

5. Try the options

6. Evaluate success

Note: If you are not successful, you may have to repeat all 6 steps, or you may not. Do I need to repeat from step 2? from step 3? From step 5?

Summary

- Troubleshooting
 - Should be performed in a linear, logical fashion
 - Know the general steps discussed
 - Keep it Simple
- These key terms discussed are important and will be used routinely in industry